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# OPERATING INSTRUCTIONS

## Venen Engel 6 Premium

Kompressionstherapie

Model: VU-IPC04

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Xiamen Weiyou Intelligent Technology Co.Ltd  
Unit 3 No. 6 Xianghong Road, Torch Hi-Tech Zone  
Industrial Park, Xiang'an District, Xiamen P.R. China.  
SRN: CN-MF-000028653

SUNGO Cert GmbH  
Harffstr. 47, 40591 Düsseldorf, Germany  
SRN: DE-AR-000010869  
(Mr.) Frank Xu  
Tel: +49 (0)211 97634133  
E-Mail: [de.rep@sunglobal.com](mailto:de.rep@sunglobal.com)

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Wellcosan GmbH  
c/o Caya Postbox 794627  
Am Börstig 5  
96052 Bamberg

Telefon: +49 (0) 666 17 48 91 05  
Mail: [service@venenengel.de](mailto:service@venenengel.de)  
Web: [www.venenengel.de](http://www.venenengel.de)

# INSTRUCTIONS FOR USE

## Venen Engel 6 Premium

Air Pressure Therapy System

Model: VU-IPC04

## PRELIMINARY REMARK

**Congratulations on the purchase of this device!** You will enjoy your Venen Engel very much! If you are satisfied with your Venen Engel, we would be very pleased if you submitted a review.

This guide includes details that should be taken into account by the user, to avoid risks, and to enable the safe use of the device.

Should you have any questions about its use, its accessories, or if you have suggestions, please do not hesitate to contact us.

**Telephone:** +49 (0) 666 17 48 91 05

**Whatsapp:** +49 1573 5990084

**Email:** [service@venenengel.de](mailto:service@venenengel.de)

**Web:** [www.venenengel.de](http://www.venenengel.de)

## BEFORE YOU GET STARTED

- Please read these instructions for use carefully before using the device.
- Keep these instructions for future consultation.
- Venen Engel may not be used if you suffer from certain health conditions. Please read the list of contraindications on the last page carefully and contact your doctor before use if you have any questions.

## PRODUCT DESCRIPTION

The compression massager consists of a device with an air pressure sensor, compressor and cuffs, which work together as one unit. The device is connected to the cuffs via a series of hoses. The compression massage is performed from the end of the limb towards the centre of the body by inflating the air chambers one after the other and then deflating them in a cycle. The pressure can be adjusted to avoid discomfort for the patient. The pressure in the cuff is controlled by a sensor and microprocessor.

In medical use, the compression massager produces sequential compression from distal to proximal, helping to improve the circulation of blood and lymph and relieve lymphoedema.

## WORKING PRINCIPLE

Compression therapy using compressed air is a therapeutic technique used in medical equipment that includes an air pump and inflatable cuffs in various designs. During use, an inflatable cuff envelops the limb to be treated. The cuff is connected to the pump via several pressure pipes. When the pump is activated, it fills the chambers of the cuff with air to apply pressure to the tissue in the limb. As the pressure decreases, the air then escapes from the chambers again. After a pause, the cycle starts again from the beginning.

The pressure on the muscles transports the venous blood to the heart and removes excess tissue fluid and metabolic products.

## SAFETY INSTRUCTIONS

The warnings and graphic symbols in this manual are intended to help you use the product safely and correctly and to avoid damage to yourself or others. The warning signs and symbols are described below:

### WARNING/CAUTION SYMBOLS

---

**Warning**

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

**Caution**

Indicates the possibility of personal injury or damage to property if the product is not used properly.

**Prohibition**

Means „prohibited“ with detailed information in words or numbers inside or next to the sign.

**Note**

Indicates the need for attention, as failure to do so may result in improper use of the product or damage to the device.

---

## EXPLANATION OF SYMBOLS USED

Store with the arrow pointing upwards

Class II equipment is equipped with double or reinforced insulation and does not require a protective conductor to be safe.

Type BF product: A medical device component that comes into contact with the patient and is electrically isolated from the patient to provide enhanced protection against leakage current.

Lot number

Serial number

Read these instructions for use carefully

Manufacturer

Date of manufacture

Keep dry

Fragile, handle with care

Do not iron

Do not wash

Do not bleach

Unique Device Identification (UDI) is a globally harmonized system for the identification and traceability of medical devices.

Keep away from direct sunlight

Range of atmospheric pressure to which the device can be safely exposed

Range of humidity to which the medical device can be safely exposed

Storage temperature range

**IP21**

Protected against solid foreign objects with a diameter of 12.5 mm (millimetres) and larger, and against vertically dripping water.

Importer

Medical device

The product complies with applicable EU directives regarding safety, health and environmental protection and has been reviewed by a notified body.

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The product is not supplied in sterile packaging

France: Recycle cardboard

Date of manufacture

France: Dispose of electronic waste at designated locations

Plastic packaging is made of polyethylene terephthalate. Please dispose of accordingly.

The cardboard is made of corrugated cardboard. Please recycle.

Made in: China

Made in: China

Authorised Representative in Switzerland

Authorized Representative in the European Union

# **INTENDED PURPOSE AND INTENDED USE**

## **INTENDED PURPOSE**

This device is intended for use by medical professionals and patients at home for the treatment of primary lymphoedema, oedema following trauma and sports injuries, oedema following immobilisation, venous insufficiencies and lymphoedema.

## **MEDICAL INDICATIONS**

Treatment of: Primary lymphoedema, oedema following trauma and sports injuries, oedema following immobilisation, venous insufficiencies and lymphoedema

## **USER**

Adults, including medical professionals and laypersons

## **PATIENT GROUP**

Adults, the patient is the intended USER

## **APPLIED PARTS**

Cuffs

## **CONTRAINDICATIONS**

- Acute pulmonary oedema
- Acute thrombophlebitis
- Acute heart failure
- Infections, deep vein thrombosis
- Episodes of pulmonary embolism wounds, lesions or tumours at or near the application site
- When increased venous and lymphatic reflux is undesirable, bone fractures or dislocations at or near the application site
- Arterial insufficiency
- Local or proximal malignancy and patients on anticoagulation therapy
- Only use during pregnancy after medical consultation

## **SIDE EFFECTS**

- Patients may experience discomfort or pain if the initial pressure is too high.
- Skin irritation may occur at the treatment site
- Fatigue or slight dizziness with prolonged use.

## SAFETY INFORMATION

### PROHIBITION

- Do not open, disassemble or modify the device under any circumstances, as this may result in fire, electric shock or other injury.
- Do not allow water or other materials (such as nails, pins and other metal objects) to enter the interior of the device.
- Position the device so that the connection to the mains can be easily disconnected.

### WARNING

- Do not use this device if you have a pacemaker, implanted defibrillator or other implanted metallic or electronic device, please consult your doctor to decide whether you are suitable for the device. Otherwise, private use may result in electrical interference or death.
- The product should not be used by infants or young children. Keep the device out of the reach of children to avoid risks such as misuse, swallowing of small parts, damage to the product or strangulation/suffocation from the power cord.
- Keep the device out of the reach of pets.
- Make sure that the plug and your hands are dry when plugging in and unplugging the power plug.
- During use of the device, you must not move around with the cuffs, but must remain seated/lying down during treatment. Keep your legs and arms still when the cuffs are on. Please empty your pockets before use, remain seated during treatment, and keep your legs and arms still while the cuffs are on.
- Use only the original accessories. The use of other accessories may cause damage to the device or failure of the treatment.
- The life of the product is 5 years. Please check that the device is fully functional before use. Using the device beyond its lifetime may result in deterioration of product performance and loss of therapeutic effect or injury. As the sensor may deteriorate, the product may no longer be able to regulate the air pressure correctly. Excessive air pressure may cause injury to the legs.
- **Emergency stop: If you experience severe pain or unusual symptoms during use, or if you want to remove the pressure in an emergency:**
  - Stop the device by pressing the on/off switch or
  - disconnect the air hose plug from the control unit or
  - remove the cuffs from your limbs.
- After storage at the maximum allowable temperature, the device should be kept at 20°C for at least 4 hours before use to ensure thermal stability.
- After storage at the minimal allowable temperature, the device should be kept at 20°C for at least 4 hours before use to ensure thermal stability.

**CAUTION** 

- Make sure that no heavy objects are placed on the air hoses, as this may cause damage to the air hose or blockage of the air flow.
- Repairs, maintenance and replacement of components must not be carried out during use. Repairs must only be carried out by specialists authorised by the manufacturer. In addition, you risk losing the warranty.
- Damage, loosening or failure of components may cause the device to malfunction. Please contact the manufacturer for repairs.
- The device should be stored in well-ventilated, dry rooms free of corrosive gases.
- Avoid shaking or dropping the unit during use or transport. Protect the unit from falls and shocks.
- Keep the device away from heat sources (e.g. radiators, cigarettes or direct sunlight) and use it only at the intended operating temperature and humidity.

## SCOPE OF DELIVERY

Control unit

Remote control

Cuffs

The scope of delivery of the cuffs depends on the set purchased.

Power cable

Air hose

A set of extension strips

## REMOTE CONTROL

**NOTE:** The remote control must be equipped with a button cell (CR2032) (size: 20 x 2.5mm DC3.0V)


## LED TOUCH SCREEN

1. Time display: LED display shows the remaining treatment time
2. Time setting: „+“ extends the treatment time; „-“ shortens the treatment time. The time range is between 1 and 99 minutes; the default setting is 30 minutes. Press „+/-“ once to add/subtract one minute. Hold „+/-“ until you reach the desired time. After the time has elapsed, the device switches off automatically.
3. Pressure display
4. Pressure setting: „+“ increases the pressure; „-“ decreases the pressure
5. Position setting: Positions 1-6 refer to the 6 chambers; the default setting includes all chambers. Press the individual chamber buttons once to turn them off (the light goes out) and the chamber will no longer operate. If you want to put the chamber back into operation, press the button again (the light comes on). While the chambers are inflating/deflating, the lights will flash.
6. On/Off button: Press once, the device is in standby mode; press again to turn off the device.
7. Start/Pause button: The original setting is on pause (green light); press the button to start the device (the light will turn blue).
8. Programme setting: There are 6 programmes available; the default setting is mode A; the other lights are off.

## PROGRAMME DESCRIPTION

- (A) Mode A: In this mode, only one chamber is inflated at a time. It starts with chamber 1 and continues to chamber 6. Then the cycle repeats.
- (B) Mode B: In this mode, chamber 1 remains inflated. Gradually, another chamber is added until all six chambers are filled with air. Then the cycle repeats.
- (C) Mode C: In this mode, chambers 1/2 are inflated simultaneously, hold the air for 2 seconds, and then start deflating; next, chambers 3/4 are inflated simultaneously, hold the air for 2 seconds, and then start deflating; then chambers 5/6 follow, hold the air for 2 seconds, and then start deflating. Then the process repeats.
- (D) Mode D: In this mode, all chambers are inflated and deflated simultaneously. Then the cycle repeats.
- (E) Mode E: Combination of B + C: sequential compression + double wave.
- (F) Mode F: Combination of A + C + D: normal + double wave + full compression.


# USE

## STEP 1: BEFORE USE

### 1. Connect air hoses

Connect the set of hoses to the cuffs, starting with the dark grey connector, and then plugging the other grey connectors onto the appropriate counterparts.

**Note:** We recommend that you do not detach the connectors afterward.

Lymph-Flow Arm


Leg cuff

### 2. Placement of the device

Ensure that the device is securely placed. Also make sure that the device is free standing and not covered by a blanket or anything similar. **Note:** Do not use the device in damp rooms or those with high humidity, such as in saunas or bathrooms.

### 3. Attaching the Multi-plug

Plug the multi-plug into the front of the device as far as it will go. Make sure to insert the multi-plug correctly into the device (the small pin into the small opening).


### 4. Combining cuffs

If you want to use accessory cuffs for simultaneous treatment of legs (leg cuff or Lymph-Flow Pants) and arms, make sure when purchasing that they have the same number of air chambers as your control unit. Only accessories with the matching number of chambers are fully compatible and functional.

If the accessories are to be operated simultaneously, use the connector provided for this purpose, which is included in the scope of delivery. Insert the connectors of both cuffs into the connector and then connect it to the opening of the main unit.

Combination air hose

## 5. Putting on the cuff

Plug the main plug of the power cord into a 230V socket. Put on the cuffs as follows:

### Arm cuffs:

An arm cuff is available for either the right or left arm. When purchasing, be sure to select the correct arm side.

**Note:** If you want to treat both arms, you will need two arm cuffs - one for the right and one for the left arm. Parallel treatment of the legs and one arm is possible. However, we recommend treating the arms one after the other and not both arms at the same time. At least one arm should always remain free during use.


### Instructions:


1. Open the clip closure of the arm cuff if it is still closed.
2. Put the arm cuff on the corresponding arm side. Pass your hand through the cuff until it rests under the armpit.
3. Pull the cuff over your shoulder.
4. Pass the closure tab behind your back under the opposite armpit and insert the end of the tab into the clip closure at the front.
5. Adjust the fit by pulling on the end of the closure tab until the cuff fits snugly and comfortably.
6. After use, loosen the arm cuff by pressing the clip closure on both sides at the same time.

### Leg cuffs and Lymph-Flow Pants:

All zippers and Velcro fasteners are open - slip your legs into the cuff.

Pull the zipper all the way up on both sides of the cuff. Tighten the Velcro fastener on both sides of the Lymph-Flow Pants - done.

### **Attaching the extension tabs (optional):**

Extensions are optionally available to increase the circumference of the leg cuffs or the Lymph-Flow Pants. When purchasing, be sure to choose the appropriate variant: For leg cuffs, you need extensions for leg cuffs, for Lymph-Flow Pants you need the corresponding extensions for Lymph-Flow Pants.

**Note:** Always make sure that the hoses are not kinked and that you are not sitting on the hoses. To avoid contamination, we recommend wearing leggings or jogging pants when using the device. Empty your pockets before use.

The zippers on the cuffs must be fully closed to prevent the cuffs from suddenly opening and causing possible damage. To prevent the zipper from opening during use, flip the slider handle towards your feet after you have pulled the pants closed. This engages a mechanical lock that prevents the Lymph-Flow pants from opening. After use, you can operate the zipper as normal to get out.

## STEP 2: CONTROL UNIT OPERATION

1. Turn the device on.
2. Set the pressure (30–240 mmHg). Start with a low setting. During the massage, you may change the pressure.
3. Time setting: Press „+“ to extend the treatment time and „-“ to shorten it. The default setting is 30 minutes. The maximum selectable time is 90 minutes.
4. Air chambers can be deactivated by pressing the corresponding button (1-6). Activated chambers are indicated by a light.
5. Select a programme. The default programme is A. The selected programme is indicated by a light.
6. Start/Pause: The original setting is the green LED display, which means the device is paused. Press the button again to make the green indicator light blue, which means the device starts working.

**Note:** The time is recalculated from zero when a time adjustment is made during operation, but an adjustment of the air pressure/position/mode and a pause during treatment do not cause the timer to restart.

### Operation using remote control

This product is equipped with a remote control. When using the remote control, the distance between the remote control and the main unit should be within 3 m. The operation of the buttons corresponds fully to that of the main unit.

**Note:** If the device is brought from a room with the minimum or maximum storage temperature to a room with a temperature of 20°C, please wait 4 hours before using the device.

### STEP 3: AFTER USE

After the set time has elapsed or after pressing the on/off switch, the device will switch off. Open the cuff. Unplug the power cord from the socket.

**Note:**

- To release the air from the air cushions more quickly, we recommend that you remove the multi-plug and press the remaining air out of the cuffs with both hands.
- It is best not to disconnect the hose from the cuffs too often. Improper use may cause damage to the pins or the hose. For storage, place the cuffs on top of each other and roll them loosely together with the attached hose set.

### STEP 4: CLEANING & STORAGE

#### 1. Cleaning

To clean the cuffs, the main device and the rest of the accessories, we recommend using a damp cloth.

Always allow the cuffs, the main unit and the accessories to dry thoroughly.

When used by different people, we recommend disinfecting the cuffs with 70% isopropanol spray.

**Note:**

- Use detergents and disinfectants free of oil, benzene, gasoline, and/or chemical agents.
- Do not wash the cuffs in the washing machine.
- Make sure that no water enters the device. If this does happen, do not use the device again until it is completely dry.
- Do not clean the device during treatment. Make sure the device is turned off before cleaning.

#### 2. Care

All cuffs are intended for compression massage and are difficult to repair. Please be careful when using.

- Do not store the device and cuff(s) near sharp or hot objects such as stoves, needles, scissors, etc.
- Store in a dry place.
- Do not store the device at low temperatures. (It may be damaged by frost).
- For long-term storage, please keep the device in the packaging.

The manufacturer has not authorised any maintenance agencies abroad. If your device has problems, please contact the distributor. The manufacturer accepts no responsibility for the results of maintenance or repair work by unauthorised persons. The user must not carry out any repairs to the device or its accessories. Please contact the distributor for repairs. Opening of the device by unauthorised parties is not permitted and will void all warranty claims. Each product has undergone systematic validation during manufacture. The performance is stable and does not require calibration or validation.

If your product does not perform as expected and the basic functions have changed during normal use, please contact the distributor. This product does not contain any parts that can be repaired or debugged by the user, so this manual does not contain technical information such as circuit diagrams and component lists. If qualified technicians of the user need them, they can request them from the manufacturer, who will provide them according to the contract.

### 3. Storage

- For long-term storage, we recommend that you store the product in the packaging you bought it in.
- Do not expose the device to direct sunlight and protect it from dirt and moisture. Store the product under the following conditions: Temperature range:  $-20^{\circ}\text{C}$  to  $55^{\circ}\text{C}$ , humidity  $\leq 93\%$ , air pressure: 50 kPa to 106 kPa.
- Store the product in frost-free places, otherwise it may be damaged. Also, be careful not to damage the cuffs during storage, especially if they are stored together with sharp objects such as scissors, or objects with sharp edges.

# **MEASURES TO PREVENT FAULTY OPERATION & EMERGENCY EVACUATION**

## **MEASURES TO PREVENT FAULTY OPERATION**

### **CAPACITIVE KEYBOARD**

Capacitive buttons are installed on the device to minimize accidental operating errors due to unintentional button presses. The touch sensitivity of these buttons is programmed in such a way that they must be pressed precisely and deliberately in order to change the pressure and / or time settings. Accidentally touching the buttons does not result in a change to the settings, nor does touching/placing your finger on them; without deliberately pressing the buttons, there is no change.

### **ACOUSTIC SIGNAL**

Every change to a setting is accompanied by a loud acoustic signal, so that you are informed acoustically of any unwanted changes and can take manual countermeasures.

### **EMERGENCY EVACUATION**

Therapy can be ended at any time by disconnecting the device from the mains. To do this, pull the plug out of the socket or press and hold the ON/OFF button until the device switches off. The cuffs deflate automatically in the event of a power failure and you can get out easily.

## TROUBLESHOOTING

### THE DEVICE DOES NOT TURN ON

If the device does not turn on, please check (A) whether the power plug has been correctly plugged into the socket and/or whether the device is turned on. If the device still does not switch on, please contact our customer service.

### THE DEVICE IS PUMPING, BUT IT ONLY BLOWS ONE OR NONE OF THE TWO CUFFS

If only one or neither of the two cuffs should inflate after starting the treatment, please check (A) whether the multi-plug is fully and correctly inserted into the device. Also, make sure (B) that the air hoses are not kinked, for example, because you are sitting on the hoses. Check (C) that the gray cuff plugs are all correctly attached to the cuff.

### THE ZIPPER OPENS DURING USE

To fix this, pull the zipper all the way up and fold the slider handle towards your feet. This engages a mechanical lock that prevents the pants from opening. After use, you can operate the zipper as you normally would to get out.

### THE CUFFS ARE INFLATING IN THE WRONG ORDER

Make sure that the multi-plug is correctly inserted into the device. The multi-plug is correctly positioned when the small nipples on the multi-plug are in the designated openings above the air outlets and the arrows on the plug are pointing upwards.

### YOU CAN HEAR AIR ESCAPING FROM THE DEVICE, HOSES, OR CUFFS


- (A) Check the hoses and plugs for damage.
- (B) Check that the multi-plug is correctly plugged into the unit.
- (C) Check if the hose is bent or has been pulled off.

## COMBINING ACCESSORIES

If you want to use other accessory cuffs, pay attention to the correct number of chambers in the cuff when buying. If you also want to use the accessory cuffs at the same time, use the connector provided for this purpose, which is included in the scope of delivery.

## USING INDIVIDUAL CUFFS

Use the dummy plug that sits in the back of the Quick Connector to block unused slots in the Quick Connector. This way you can also use just a single cuff (e.g. an arm cuff) without having to connect a second cuff to the quick connector.

## DISPOSAL

At the end of the product life cycle, do not dispose of this product in normal household waste, but take it to a collection point for recycling electronic equipment. Waste electrical and electronic equipment can potentially have harmful effects on the environment. Improper disposal can cause harmful toxins to accumulate in the air, water, and soil, affecting human health. You are responsible for the proper disposal of the device. For information on disposal, please contact your local authority or your distributor.

## TECHNICAL DETAILS

<b>Number of chambers</b>	6
<b>Size</b>	30 x 23,7 x 12,6 cm
<b>Weight (with all accessories)</b>	Approx. >5 kg
<b>Time span</b>	1-90 minutes
<b>Pressure range</b>	30-240 mmHg
<b>Programme</b>	6
<b>Application site</b>	Legs, arms
<b>Input voltage</b>	220-230V, 50-60Hz
<b>Power consumption</b>	65W
<b>Volume</b>	≤65dBa (Margin of error: ± 3dBa)
<b>Type of protection against electric shock</b>	Class II
<b>Operating mode</b>	Continuous operation
<b>Waterproof</b>	IP21
<b>Product lifetime</b>	5 years
<b>Operating environment</b>	10-40°C, 30-85%RH, 70-106kPa
<b>Storage environment</b>	-40-70°C, 10-100%RH, 50-106kPa
<b>Material composition (cuff)</b>	50% nylon, 50% TPU (polyurethane)

## REMOTE CONTROL

<b>Power supply</b>	DC 3.0V, 1x CR2032 cell
<b>Transmission type</b>	Infrared
<b>Range</b>	$\leq 3\text{m}$
<b>Size</b>	92,8 x 48,2 mm
<b>Weight</b>	20g

## IMPORTANT INFORMATION REGARDING ELECTROMAGNETIC COMPATIBILITY (EMC)

The list of cables and their maximum length is as follows:

Cable name	Cable length	Weather protection
Power cable	1.9 m	No

### WARNING

The use of mobile phones or microwave ovens, RF surgical equipment, magnetic resonance imaging equipment or other radio radiation equipment near this product may cause malfunction or loss of essential performance characteristics, impairing measurement accuracy.

The use of this device in close proximity to or stacked with other devices should be avoided as this may cause malfunction. If such use is necessary, this device and the other devices should be observed to ensure that they are operating properly.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) to any part of the air pressure therapy system, including the cables specified by the manufacturer. Otherwise, the performance of the device may be degraded.

The use of accessories, transducers and cables not specified or provided by the manufacturer of this device may result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and cause improper operation.

**CAUTION** 

Security, anti-theft and radio frequency identification (RFID) devices. Some electromagnetic anti-theft systems and metal detectors, such as those used at entrances and exits of department stores, libraries and other public places, as well as security checks at airports, can interfere with the compression device. In addition, RFID devices, which are often used to read ID cards, as well as some label deactivation devices, such as those used at checkouts in shops and loan counters in libraries, can also interfere with the compression device.

Please do not use the device near these locations. If you must pass one of these devices, turn off your compression device. Check the status of your compression device before each use to ensure that it is functioning normally.

The use of shortwave diathermy, microwave diathermy, or therapeutic ultrasound diathermy and electrocautery devices near this product may cause malfunctions or loss of essential performance characteristics.

<b>Guidelines and manufacturer's declaration – electromagnetic emissions</b>		
The device is intended for use in the electromagnetic environment specified below. The customer of the device user should ensure that it is used in such an environment.		
<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment - Guidance</b>
Conducted and radiated RF emissions CISPR 11	Group 1 Class	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted RF emissions CISPR 11	Group 1 Class B	The compression therapy device is suitable for use in all facilities, including households and those directly connected to the public low-voltage power supply that supplies the building for residential purposes, with the exception of proximity to active RF surgical equipment and the RF-shielded room for magnetic resonance imaging.
Radiated RF emissions CISPR 11	Group 1 Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

## Guidelines and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below.  
The customer or user of the device should ensure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - Guide
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact; ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact; ±2kV, ±4kV, ±8Kv, ±15 kV air	Floors should be made of wood, concrete, or ceramic tiles. If floors are covered with synthetic material, relative humidity should be at least 30%.
Radiated RF EM fields IEC 61000-4-3	3 V/m (professional health-care facilities); 10 V/m (domestic care), 80 MHz-2.7 GHz 80% AM at 1 kHz	10V/m (professional healthcare facility) 80 MHz-2.7 GHz 80% AM at 1 kHz	
Electrical fast transients/ bursts IEC 61000-4-4	±2 kV AC for power supply lines; ±1 kV DC power/signal lines. 100 kHz repetition frequency	±2kV for power supply lines (AC)	The mains power quality should correspond to that of a typical commercial or clinical environment.
Waveform IEC 61000-4-5	±0.5 kV, ±1 kV line(s) to line(s); ±0.5 kV, ±1 kV line(s) to ground	±0.5kV, ±1kV Line(s) to line(s)	The mains power quality should be typical of a commercial or clinical environment.
Line-borne disturbances induced by RF fields IEC 61000-4-6	3V 0.15MHz-80MHz, 6V in ISM bands between 0.15 MHz and 80 MHz (professional healthcare facilities), 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz	Applicable	
<p><b>Note:</b> The ISM (industrial, scientific, and medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz, and 50.0 MHz to 54.0 MHz.</p>			
Radiated magnetic fields with mains frequency IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz	Magnetic fields at power frequency should be at levels characteristic of a typical location in a typical commercial or clinical environment.

Voltage dips IEC 61000-4-11	0% $U_T$ , 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°; 0% $U_T$ , 1 cycle and 70% $U_T$ , 25/30 cycles Single-phase: at 0°	Applicable	The mains power quality should be typical of a commercial or clinical environment.
Voltage interruptions IEC 61000-4-11	0% $U_T$ , 250/300 cycles	Applicable	
 <b>NOTE:</b> $U_T$ is the AC mains voltage before the test voltage is applied. Example: 25/30 means 25 periods at 50 Hz or 30 periods at 60 Hz.			

#### Guidelines and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below.  
The customer or user of the device should ensure that it is used in such an environment.

Immunity test	IEC 60601-1-2 Test level	Compliance level	Electromagnetic environment - Guide
Near-field areas of RF radio communication equipment IEC 61000-4-3	See the following table	Corresponds	

**Test specifications for enclosure port immunity to RF wireless communications equipment**

Test frequency (MHz)	Band 1 (MHz)	Service ·1	Modulation <sup>b)</sup>	Maximum power (W)	Distance (m)	Immunity Test Level (V/m)
385	380-390	TETRA 400	Pulse modulation <sup>b)</sup> 18 Hz	1,8	0,3	27
450	430-470	GMRS460, FRS 460	FM <sup>a)</sup> ±5 kHz deviation 1kHz sine	2	0.3	28
710 745 780	704-787	LTE Band 13, 17	Pulse modulation <sup>b)</sup> 217 Hz	0.2	0.3	9
810 870 930	800-960	GSM800/900, TETRA 800, iDEN 820, CDMA850, LTE Band 5	Pulse modulation <sup>b)</sup> 18 Hz	2	0.3	28
1720 1845	1700-1990	GSM1800; CDMA 1900;	Pulse modulation <sup>b)</sup> 217 Hz	2	0.3	28
1970		GSM 1900; DECT; LTE Band 1,3, 4, 25; UMTS	Pulse modulation <sup>b)</sup> 217 Hz			
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation <sup>b)</sup> 217 Hz	2	0.3	28
5240 5500 5785	5100-5800	WLAN 802.11 a/n	Pulse modulation <sup>b)</sup> 217 Hz	0.2	0.3	9

**NOTE:** If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME DEVICE or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted according to IEC 61000-4-3.

- a. For some services, only the uplink frequencies are included.
- b. The carrier shall be modulated using a 50% duty cycle square wave signal.
- c. As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be the worst case.

The air pressure therapy system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the equipment can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the air pressure therapy system.

**Recommended separation distances between portable and mobile RF communications equipment and air pressure therapy equipment**

Rated maximum output power of transmitter/W.	Separation distance according to frequency of transmitter/m		
	150kHz ~ 80MHz $d = 1.2 \sqrt{P}$	80MHz ~ 800MHz $d = 1.2 \sqrt{P}$	800MHz ~ 2.5GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1:** At 80 MHz, the separation distance for the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## WARRANTY

We offer a 24-month warranty on the control unit. In the event of malfunction, it may be necessary to send the device back for inspection. Please ensure that the shipping box is padded to avoid possible damage during transport. Unfortunately, no claims under the warranty can be brought for defects caused by improper shipping.

**Terms of warranty:** insofar as a legal warranty obligation exists, a replacement will either be provided free of charge or the manufacturer shall have the option to repair the device, under exclusion of the right of cancellation or a reduction in the purchase price. If repeated attempts to effect repair are unsuccessful, or a replacement unit also proves to be defective due to the manufacturer's fault, the customer shall be entitled to cancel the purchase or receive a reduction in the purchase price.

**Exclusion of warranty:** Excluded from the warranty are defects due to damage caused by the use of force, improper operation, externally applied force, or modifications and repair work performed by third parties, e.g. defects caused by incorrectly rated or short-circuited fuses, or defects attributable to normal wear and tear.

**Conditions of warranty:** The warranty claim exists only in connection with the original proof of purchase (invoice). Please keep the proof of purchase in a safe place. Should you have any problems or questions when using the system, please do not hesitate to contact us at +49 (0) 666 174 891 05 or by mail at [service@venenengel.de](mailto:service@venenengel.de).

## OTHER INFORMATION

Xiamen Weiyou Intelligent Technology Co.Ltd  
Unit 3 No. 6 Xianghong Road, Torch Hi-Tech Zone  
Industrial Park, Xiang'an District, Xiamen P.R. China.  
SRN: CN-MF-000028653

SUNGO Cert GmbH  
Harffstr. 47, 40591 Düsseldorf, Germany  
SRN: DE-AR-000010869  
(Mr.) Frank Xu  
Tel: +49 (0)211 97634133  
Email: [de.rep@sunglobal.com](mailto:de.rep@sunglobal.com)

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